Química Analítica I

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Material Safety Data Sheet

Version 4.1
Revision Date 10/23/2010
Print Date 02/08/2011

1. PRODUCT AND COMPANY IDENTIFICATION

Product name	Chromium(III) acetate hydroxide		
Product Number	318108		
Brand	Aldrich		
Product Use	For laboratory research purposes.		
Supplier	Sigma-Aldrich Canada, Ltd. 2149 Winston Park Drive OAKVILLE ON L6H 6J8 CANADA	Manufacturer	Sigma-Aldrich Corporation 3050 Spruce St. St. Louis, Missouri 63103 USA
Telephone	+19058299500		
Fax	+19058299292		
Emergency Phone # (For both supplier and manufacturer)	1-800-424-9900		
Preparation Information	Sigma-Aldrich Corporation Product Safety - Americas Region 1-800-521-8956		

2. HAZARDS IDENTIFICATION

Emergency Overview

WHMIS Classification

Not WHMIS controlled.

Not WHMIS controlled.

GHS Classification

Acute toxicity, Inhalation (Category 4)

Acute toxicity, Dermal (Category 4)

Acute toxicity, Oral (Category 4)

Skin irritation (Category 2)

Eye irritation (Category 2A)

Specific target organ toxicity - single exposure (Category 3)

GHS Label elements, including precautionary statements

Pictogram

Signal word

Warning

Hazard statement(s)

H302 + H312

Harmful if swallowed or in contact with skin.

H315

Causes skin irritation.

H319

Causes serious eye irritation.

H332

Harmful if inhaled.

H335

May cause respiratory irritation.

Precautionary statement(s)

P261

Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.

P264

Wash skin thoroughly after handling.

P270

Do not eat, drink or smoke when using this product.

P271

Use only outdoors or in a well-ventilated area.

P280

Wear protective gloves/ eye protection/ face protection.

P301 + P312

IF SWALLOWED: Call a POISON CENTER or doctor/ physician if you feel unwell.

Aldrich - 318108

Page 1 of 6

(Material Safety Data Sheet, MSDS o Safety Data Sheet, SDS)

1. PRODUCT AND COMPANY IDENTIFICATION

Product name	Chromium(III) acetate hydroxide		
Product Number	318108		
Brand	Aldrich		
Product Use	For laboratory research purposes.		
Supplier	Sigma-Aldrich Canada, Ltd 2149 Winston Park Drive OAKVILLE ON L6H 6J8 CANADA	Manufacturer	Sigma-Aldrich Corporation 3050 Spruce St. St. Louis, Missouri 63103 USA
Telephone	+19058299500		
Fax	+19058299292		
Emergency Phone # (For both supplier and manufacturer)	1-800-424-9300		
Preparation Information	Sigma-Aldrich Corporation Product Safety - Americas Region 1-800-521-8956		

2. HAZARDS IDENTIFICATION

Emergency Overview

WHMIS Classification

Not WHMIS controlled.

Not WHMIS controlled.

GHS Classification

Acute toxicity, Inhalation (Category 4)

Acute toxicity, Dermal (Category 4)

Acute toxicity, Oral (Category 4)

Skin irritation (Category 2)

Eye irritation (Category 2A)

Specific target organ toxicity - single exposure (Category 3)

GHS Label elements, including precautionary statements

Pictogram

Signal word

Warning

Hazard statement(s)

H302 + H312

Harmful if swallowed or in contact with skin.

H315

Causes skin irritation.

H319

Causes serious eye irritation.

H332

Harmful if inhaled.

H335

May cause respiratory irritation.

Precautionary statement(s)

P261

Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.

P264

Wash skin thoroughly after handling.

P270

Do not eat, drink or smoke when using this product.

P271

Use only outdoors or in a well-ventilated area.

P280

Wear protective gloves/ eye protection/ face protection.

P301 + P312

IF SWALLOWED: Call a POISON CENTER or doctor/ physician if you feel unwell.

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
Dr. Dimbrij Rylvén
Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.

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GHS PICTOGRAMS

Health Hazard Carcinogens, respiratory sensitisers, reproductive toxicity, target organ toxicity, germ cell mutagens 	Flame Flammable gases, liquids, & solids; self-reactives; pyrophorics; 	Exclamation Mark Irritant, dermal sensitiser, acute toxicity (harmful) 
Gas Cylinder Compressed gases; liquefied gases; dissolved gases 	Corrosion Skin corrosion; serious eye damage 	Exploding Bomb Explosives, self-reactives, organic peroxides 
Flame Over Circle Oxidisers gases, liquids and solids 	Environment Aquatic toxicity 	Skull & Crossbones Acute toxicity (severe) 

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The 12 Principles of
GREEN CHEMISTRY

Green chemistry is an approach to chemistry that aims to maximize efficiency and minimize hazardous effects on human health and the environment. While no reaction can be perfectly 'green', the overall negative impact of chemistry research and the chemical industry can be reduced by implementing the 12 Principles of Green Chemistry wherever possible.

1. WASTE PREVENTION  Prioritize the prevention of waste, rather than cleaning up and treating waste after it has been created. Plan ahead to minimize waste at every step.	7. USE OF RENEWABLE FEEDSTOCKS  Use chemicals which are made from renewable (i.e. plant-based) sources, rather than other, equivalent chemicals originating from petrochemical sources.
2. ATOM ECONOMY  Reduce waste at the molecular level by maximizing the number of atoms from all reagents that are incorporated into the final product. Use atom economy to evaluate reaction efficiency.	8. REDUCE DERIVATIVES  Minimize the use of temporary derivatives such as protecting groups. Avoid derivatives to reduce reaction steps, resources required, and waste created.
3. LESS HAZARDOUS CHEMICAL SYNTHESIS  Design chemical reactions and synthetic routes to be as safe as possible. Consider the hazards of all substances handled during the reaction, including waste.	9. CATALYSIS  Use catalytic instead of stoichiometric reagents in reactions. Choose catalysts to help increase selectivity, minimize waste, and reduce reaction times and energy demands.
4. DESIGNING SAFER CHEMICALS  Minimize toxicity directly by molecular design. Predict and evaluate aspects such as physical properties, toxicity, and environmental fate throughout the design process.	10. DESIGN FOR DEGRADATION  Design chemicals that degrade and can be discarded easily. Ensure that both chemicals and their degradation products are not toxic, bioaccumulative, or environmentally persistent.
5. SAFER SOLVENTS & AUXILIARIES  Choose the safest solvent available for any given step. Minimize the total amount of solvents and auxiliary substances used, as these make up a large percentage of the total waste created.	11. REAL-TIME POLLUTION PREVENTION  Monitor chemical reactions in real-time as they occur to prevent the formation and release of any potentially hazardous and polluting substances.
6. DESIGN FOR ENERGY EFFICIENCY  Choose the least energy-intensive chemical route. Avoid heating and cooling, as well as pressurized and vacuum conditions (i.e. ambient temperature & pressure are optimal).	12. SAFER CHEMISTRY FOR ACCIDENT PREVENTION  Choose and develop chemical procedures that are safer and inherently minimize the risk of accidents. Know the possible risks and assess them beforehand.

  © COMPOUND INTEREST 2015; WWW.COMPOUNDCHEM.COM
Shared under a CC Attribution-NonCommercial-NoDerivatives licence.

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
The diagram shows a sequence of four steps for correcting a mistake in a notebook:

- Step 1.** Draw a single line through the error. (Illustrated by a line through the number 28 in a list: 25, 26, ~~28~~, 27, 28).
- Step 2.** Make the correction next to the error with an asterisk and a number. (Illustrated by writing '27 *1' next to the crossed-out 28).
- Step 3.** Write an explanation for the error making reference to the asterisk number. (Illustrated by writing '*1 wrong entry Type error' below the list).
- Step 4.** Sign and date the correction. (Illustrated by writing 'R. Cayueh' and '11-27-20' at the bottom right).

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